

You can make a difference!

GC Tooth Mousse contains 10% (wt) RECALDENT™ CPP-ACP to provide you with a concentrated source of bio-available calcium and phosphate for strengthening and protection of teeth. And it tastes delicious!

After brushing, treat your teeth with GC Tooth Mousse and help everyone in the family achieve a healthy smile.




What is RECALDENT™ CPP-ACP?


Nature has developed highly efficient protein systems in saliva, blood and milk that stabilise and transport calcium and phosphate which are the fundamental building blocks of our teeth and bones. Cow's milk is recognised as having the most effective protein system for delivery of calcium and phosphate in a bio-available form, and extensive research over the past 25 years has identified and refined the specific peptide found in milk which so efficiently transports these essential minerals to help build, strengthen and protect our teeth – casein phosphopeptide amorphous calcium phosphate – otherwise known as RECALDENT™ (CPP-ACP).



So easy to use!

- 

1 Apply a small amount of GC Tooth Mousse to all teeth with your finger.
- 

2 Use your tongue to spread the mousse evenly around your teeth - including hard to reach
- 

3 For maximum benefit, leave GC Tooth Mousse on the teeth for as long as possible. The minimum recommended application time is three minutes. At the end of the application you can expectorate, rinse or swallow the remainder.

GC Tooth Mousse comes in individual 40g tubes in a variety of delicious flavours.



GC Asia Dental Pte. Ltd
Changi Logistics Centre
19 Loyang Way #06-27 Singapore 508724
Tel + 65 6546 7588 Fax + 65 6546 7577
www.gcasia.info © GC Corporation 2010



GC Tooth Mousse for Children



Nature knows best

Made from milk

GC Tooth Mousse contains RECALDENT™ (CPP-ACP), a unique ingredient developed at The School of Dental Science, The University of Melbourne, Victoria, Australia. RECALDENT and RECALDENT Device are trademarks used under licence. GC Tooth Mousse should not be used by people with milk protein allergies. If any allergic reaction occurs, this may indicate sensitivity to the benzoate preservatives, or to some other component of the product. In this event, discontinue use of the product and contact your physician.



Baby Mousse

GC Tooth Mousse is especially useful for under 2's where it leaves a surface film that raises levels of essential minerals and well as inhibiting the growth of decay causing bacteria.

“Newly erupted teeth have yet to complete their enamel maturation and until this occurs they are more vulnerable to acid attack. Boosting levels of calcium and phosphate in the saliva facilitates the normal post eruption maturation process and replaces mineral loss on a daily basis. RECALDENT™ (CPP-ACP) is derived from cows’ milk and is ideal for protecting deciduous teeth at a time when oral care is difficult.”

PROF LJ WALSH, UNIVERSITY OF QUEENSLAND



Junior Mousse

This is the ideal time to create good oral care habits which will have benefits for a lifetime. It is also a time when permanent teeth are erupting and additional protection, especially of erupting molars may be required.

“Starting early means healthy teeth and a healthy lifestyle for the future.”

PROF LJ WALSH, UNIVERSITY OF QUEENSLAND

Mousse for the mum to be



“During pregnancy, hormonal changes result in a dramatic decrease in the ability of saliva to counteract acid attack. At the same time, nausea from morning sickness can lead to frequent reflux of acid contents of the stomach, particularly in the second and third months of pregnancy. Reduced protection from saliva, combined with frequent acid reflux, means that many pregnant women are at high risk of tooth decay during pregnancy. Regular use of GC Tooth Mousse during pregnancy helps to maintain elevated levels of calcium and phosphate in the saliva, protecting the teeth from dental erosion. GC Tooth Mousse will also help soothe sensitive teeth and protect areas of decay. There is also an additional benefit because GC Tooth Mousse will inhibit growth of caries-causing bacteria which can be easily transmitted to the newborn child.”

PROF LJ WALSH, UNIVERSITY OF QUEENSLAND

Ortho Mousse



“GC Tooth Mousse can help prevent plaque accumulation around brackets, archwires, springs and other appliances, and can assist the saliva in countering acids produced by dental plaque. A regular treatment of GC Tooth Mousse throughout the full course of orthodontic treatment can prevent areas of decalcification (white spots) developing.”

PROF LJ WALSH, UNIVERSITY OF QUEENSLAND

Application guide

Indication	How often	Duration
General prevention	Once in the evening after flossing and brushing.	Ongoing as part of a general prevention program
Additional comments: Regular use of GC Tooth Mousse will help maintain a healthy oral environment.		
Active decay	Twice daily after flossing and brushing.	Until risk of future decay has been reduced.
Additional comments: Your dental professional can assess and determine the potential source(s) of risk and can suggest how to reduce risk.		
White spot lesions	Twice daily after flossing and brushing.	8-12 weeks and thereafter required.
Additional comments: There are different types of white spots and your dental professional may need to pre-treat the white spot surface prior to starting application of GC Tooth Mousse.		
Orthodontic treatment	Twice daily after flossing and brushing.	During entire orthodontic process.
Additional comments: Complete treatment 12 weeks after the finish of orthodontic procedures or thereafter as required to reverse any white spot lesions.		
Molar Incisor Hypo-mineralisation (MIH)	Minimum twice a day, after flossing and brushing, and as required for sensitivity.	Continuous treatment as required.
Additional comments: Depending on the severity of MIH your dental professional may need to undertake additional procedures to protect these teeth.		